



Ionenanalytik • Ion analysis • Analyse des ions • Análisis iónico

Certificate

Metrohm AG CH-9101 Herisau
Switzerland

This is to certify that the instrument below conforms to the specifications mentioned.

Instrument number: **3.496.8800**
Instrument designation: PT100 PRECISION RESISTOR SET
Serial number: 08147
Inventory number: -
Calibration Master: FLUKE 8508A PM10/643
Serial number: 913152733
Certified by: SCS Cal.Serv.Reg.No.02
Certificate number: 12534

The Resistor Set includes the following resistors:

Type	Measured Value [Ω]
3496881008154	119.9931
3496882008158	149.9871
3496883008160	199.9947
3496884007152	119.9829
3496885007134	200.0094

Values measured at a room temperature of $25.3\text{ }^{\circ}\text{C} \pm 2\text{ }^{\circ}\text{C}$

Technical data:
Temperature coefficient: 2ppm/ $^{\circ}\text{C}$
Longterm stability: 0.025‰/year

Date:

Tested by: R. Helfenberger

19.09.2012

Next calibration: